

# TYMPANIC THERMOMETER TESTER

*The TTT is a new innovation in biomedical test equipment. It has been designed to assure the correct performance of tympanic thermometers.*

*This new product is of great benefit to 'In-house' Bio-med's, Third Party Service Providers, Hospital Maintenance Dept's. & O.E.M's, as it brings many advantages, for example improving productivity, saving money and of course, compliance with AS:3551:2002 - In-Service Testing of Essential Performance Parameters.*

*Up until now it was difficult to confirm accuracy, as a multi-mode, portable thermometer tester was not available. Now the TTT will eliminate doubt and provide peace of mind. Clinical staff can have more confidence in Tympanic Thermometer readings, as suspect units can be tested in seconds with high accuracies.*

*The TTT is operated on battery or mains power, hence is very portable and efficient to use. Tympanic thermometers can be tested on site, eliminating unnecessary down-time.*

*With seven different operating modes, including modes for the most popular tympanic thermometers available, the TTT is gives the user a great deal of flexibility.*

*It incorporates three, infrared radiation sources, to test tympanic thermometers at different temperatures: 34.0, 37.0 & 40.0 degrees Celsius. The TTT can verify if the tympanic thermometer under test, is functioning correctly, with an accuracy of +/- 0.1 Degrees. Calibration is traceable to national Australian standards.*

*...How accurate are your thermometers ???*



Supplied by: